

Lithotripsy Baskets and Spirals

Single Use



LEF-GML-SU-00-EN / LL-0385 Vers. 1.0 (11/2020)

Combination Concept

Components for planned and emergency lithotripsies

Simple and Safe Handling

Convenient handling for the assembly and detachability of components

Pre-assembled Set Configurations

For specific treatment requirements



ENDO-FLEX GmbH

Alte Hünxer Str. 115 • 46562 Voerde • Germany
Tel.: +49 281 94400-0 • Fax: +49 281 94400-11
Email: sales@endo-flex.de • www.endo-flex.de

Lithotripsy Baskets – Single Use


Lithotripsy Baskets – Single Use



Lithotripsy Basket – 4 Braided Wire

Suitable for lithotripsy handle type I & III; tube length: 195 cm; total system length: 400 cm

| Catalogue No. | Description | Sheath Outer Diameter (mm) | Suitable for Lithotripsy Spiral Ø (mm) | Basket Width in Opened Condition (mm) | Basket Length in Stretched Condition (mm) |
|---------------|--------------------|----------------------------|----------------------------------------|---------------------------------------|-------------------------------------------|
| E122264N | 4-Wire Basket Type | 2.6 | 2.6 | ~20 | 40 |
| E122265N | 4-Wire Basket Type | 2.6 | 2.6 | ~25 | 50 |
| E122266N | 4-Wire Basket Type | 2.6 | 2.6 | ~30 | 60 |
| E122267N | 4-Wire Basket Type | 2.6 | 2.6 | ~35 | 70 |
| E122336N | 4-Wire Basket Type | 2.9 | 3.3 | ~30 | 60 |





Lithotripsy Basket – 6 Braided Wire

Suitable for lithotripsy handle type I & III; tube length: 195 cm; total system length: 400 cm

| Catalogue No. | Description | Sheath Outer Diameter (mm) | Suitable for Lithotripsy Spiral Ø (mm) | Basket Width in Opened Condition (mm) | Basket Length in Stretched Condition (mm) |
|---------------|--------------------|----------------------------|----------------------------------------|---------------------------------------|-------------------------------------------|
| E124263N | 6-Wire Dormia Type | 2.6 | 2.6 | ~12 | 30 |
| E124264N | 6-Wire Dormia Type | 2.6 | 2.6 | ~18 | 40 |
| E124265N | 6-Wire Dormia Type | 2.6 | 2.6 | ~25 | 50 |
| E124266N | 6-Wire Dormia Type | 2.6 | 2.6 | ~30 | 60 |
| E124267N | 6-Wire Dormia Type | 2.6 | 2.6 | ~35 | 70 |
| E124335N | 6 Wire Dormia Type | 2.9 | 3.3 | ~25 | 50 |
| E124337N | 6-Wire Dormia Type | 2.9 | 3.3 | ~35 | 70 |
| E128336N | 6-Wire Basket Type | 2.9 | 3.3 | ~30 | 60 |



Fitting Lithotripsy Handles – Mechanical Lithotripsy

| Catalogue No. | Description |
|---------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| 120000N | <p>Lithotripsy Handle Type I</p> <p>Reusable, with coaxial fixation, ergonomic hand wheel suitable for planned Lithotripsy only</p>  |
| 120600 | <p>Lithotripsy Handle Type III</p> <p>Reusable, user-friendly twist & pull technology, 1-push safety system suitable for planned and emergency Lithotripsy</p>  |

Lithotripsy Basket - 4 Braided Wire

Suitable for lithotripsy handle type I; compatible with guide wire .035 inch; tube length: 195 cm; total length: 400 cm



| Catalogue No. | Description | Sheath Outer Diameter (mm) | Suitable for Lithotripsy Spiral Ø (mm) | Basket Width in Opened Condition (mm) | Basket Length in Stretched Condition (mm) |
|---------------|--------------------|----------------------------|----------------------------------------|---------------------------------------|-------------------------------------------|
| E122263N-GW | 4-Wire Basket Type | 2.6 | 2.6 | ~ 15 | 30 |
| E122265N-GW | 4-Wire Basket Type | 2.6 | 2.6 | ~ 25 | 50 |
| E122266N-GW | 4-Wire Basket Type | 2.6 | 2.6 | ~ 30 | 60 |
| E122267N-GW | 4-Wire Basket Type | 2.6 | 2.6 | ~ 35 | 70 |
| E122268N-GW | 4-Wire Basket Type | 2.6 | 2.6 | ~ 40 | 80 |



Lithotripsy Basket – 6 Braided Wire

Suitable for lithotripsy handle type I; compatible with guide wire .035 inch; tube length: 195 cm; total length: 400 cm



| Catalogue No. | Description | Sheath Outer Diameter (mm) | Suitable for Lithotripsy Spiral Ø (mm) | Basket Width in Opened Condition (mm) | Basket Length in Stretched Condition (mm) |
|---------------|--------------------|----------------------------|----------------------------------------|---------------------------------------|-------------------------------------------|
| E124266N-GW | 6-Wire Dormia Type | 2.6 | 2.6 | ~ 30 | 60 |



Lithotripsy Spirals – Single Use

Lithotripsy Spirals – Single Use



Lithotripsy Spiral

Made of stainless steel

| Catalogue No. | Suitable for Lithotripsy Handle | Outer Diameter (mm) | Total Length (cm) |
|---------------|---------------------------------|---------------------|-------------------|
| E120261 | Type I | 2.6 | 195 |
| E120602 | Type III | 2.6 | 195 |



Emergency Lithotripsy Spiral

Made of stainless steel

| Catalogue No. | Suitable for Lithotripsy Handle | Outer Diameter (mm) | Total Length (cm) |
|---------------|---------------------------------|---------------------|-------------------|
| E120604 | Type III | 4.0 | 100 |





Lithotripsy Set Type III

Consists of storage case, handle, spirals and baskets

| Catalogue No. | Suitable for Lithotripsy Handle | Outer Diameter (mm) | Component Catalogue No. |
|---------------|---------------------------------|----------------------------------|-------------------------|
| E120610 | Type III | 1 x Lithotripsy Handle Type III | 120600 |
| | | 1 x Lithotripsy Storage Case | 120601 |
| | | 3 x Lithotripsy Spirals | E120602 |
| | | 1 x Emergency Lithotripsy Spiral | E120604 |
| | | 1 x Lithotripsy Basket | E120265N |
| | | 1 x Lithotripsy Basket | E120266N |
| | | 1 x Lithotripsy Basket | E120267N |

Fitting Lithotripsy Handles – Mechanical Lithotripsy

| Catalogue No. | Description |
|---------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| 120000N | <p>Lithotripsy Handle Type I</p> <p>Reusable, with coaxial fixation, ergonomic hand wheel suitable for planned Lithotripsy only</p>  |
| 120600 | <p>Lithotripsy Handle Type III</p> <p>Reusable, user-friendly twist & pull technology, 1-push safety system suitable for planned and emergency Lithotripsy</p>  |